


Cassette Plasma Sterilizer



Cartridge




Printer



Chamber



Incubator




Foot sensor switch

Model	Capacity	Chamber size (cm)	Machine size (cm)	Power	Voltage
SQ-D-220	190 L	82×51×46	112×83×160	5.0 KW	220V 50HZ

Consumables



Cartridge



Tyvek Pouch



Non-woven Sheet



Biological Indicator



Indicator tape



Chemical Strip



H2O2 Capsule



Printer paper

Performance

Power supply	AC 220V 50 Hz or AC 110V 60 Hz
Language	English / Chinese or Spanish/Chinese
Cartridge	12 pcs capsule
H2O2 consumption	1.Short cycle 1 capsule . 2. Long cycle 2 capsule . 3. H2O2 solution valid time 1 year .
Lifting Door	Mandril-driving auto overhead door with reliable sealing function and safe lock device to guarantee personal safety.
sterilization temperature	45-55°C
Control system	Siemens PLC , 7" HMI touch panel The process is fully automatic
Record	Built-in thermal printer , real-time print out data USB port available for historical data download
Foot sensor switch	Open or shut down the chamber door by foot , when operator hands occupied .

Exclusive Distributor



PT. PERMANA PUTRA MANDIRI
 MEDICAL EQUIPMENTS SUPPLIES AND SERVICES
 Jl. Taman Mini Indonesia Indah (TMII)
 Pintu II Atas No. 43 Cipayung, Jakarta Timur, Indonesia
 Telp. : +62-21 841 5976, Fax: +62-21 840 7865
 Email : info@ptppm.id, Website : www.ptppm.id