

High Frequency Inverter X-ray System



PT. Permana Putra Mandiri
Medical Equipment Supplies and Services
Jl. Taman Mini Indonesia Indah (TMII) Pintu II atas No. 43 Cipayung, Jakarta Timur - Indonesia
Phone : 021 841 5976, Fax: 021 840 7865
Email: info@ptppm.id, marketing@alatkedokteran.id
www.ptppm.id | www.alatkedokteran.id





Can reduce shooting time about 50% than general method when

WSR-30 system films human body by high frequency inverter method and because control console is small and light, is possible by wall tapestry and Radiographic is available by remote control switch conveniently that user films. And have A.P.R function to basis.

SPECIFICATIONS

Component	Item	WSR-30
Consol & Generator	KVP Range	40~125KV
	mA Range	50mA~300mA
	mAs Range	0.5mAs~500mAs
	Time Range	0.01~5sec
	Sys-Self	LED code display
	Diagnosis	
	A.P.R	576kinds
	Output Frequency	60KHZ
Tube unit	Tube Voltage	500mA/125KVP
	Anode Heat	140,000HU
Tube Stand	U/D mov. dis	Approx. 1300mm(680 to 1980mm)
	L/D mov. dis	Approx. 1480mm Rail 2400mm
	Revo. angle	360°
Table	L/R mov. dis	±350mm (2,000L×728W×710H)
Collimator	Lighting	More than 160 Lux at SID 100cm
H.T cable	KVP Range	Max. 150KVP
Option	Bucky Stand	
	Ceiling	
	6way-table	

