
Feature

- 12 leads synchronous acquisition, display on the same screen, synchronization print
- 9 inch color LCD screen, support touch screen operation
- Manual/Auto printing: 3CH+3R, 6CH, 6CH+1R, 12CH, 5 minutes arrhythmia analysis of waveform
- 60 seconds rhythm analysis of waveform
- 200 ECG data storage, SD card supported
- Automatic analysis, calculation, waveform review
- Alarm of electrodes breaking off, Human diagram indicates the unconnected lead
- V1,V3,V5 chest lead examination mode, fit for massive physical examination
- Built-in rechargeable lithium battery, AC/DC power conversion
- Support ECG-20NT digital ECG work station, achieve depth analysis of ECG



Strong and lightweight trolley





 WUHAN ZONCARE BIO-MEDICAL ELECTRONICS CO., LTD.

 Add. Zoncare Building, No. 390, High-tech 2nd Road, Eastlake High-tech Development Zone, Wuhan, Hubel 430206 P. R. China.

 Tei: +88-27-8777 0581

 Tei: +88-27-8777 0581

 Feizi: Info@zoncare.com

 http://www.zoncare.com

 2015 Wuhan Zoncare Bio-Medical Electronics Co. Ltd. Al Lights reserved. Specifications subject to change

http://www.zoncare.com





